

• *biol*atex

“the *latex*
reagent company”



⋮ *biollatex*

WWW.BIOLATEX.COM



• *biol*atex develops and manufactures - OEM, bulk and packed - liquid reagents ready to use for latex enhanced immunoturbidimetric methods.

The covalent bond between the antigens or antibodies and the micro-particles of latex provides our reagents an excellent accuracy, a high reliability and an up to two years stability.

Strong performance is also assured through the time, following strict manufacturing quality guidelines that allow our reagents obtain similar responses irrespectively of the production batch.

Our commitment with the lab instrument companies and distributors is to meet their needs and establish, if necessary, special contracts for the development and manufacture of reagents.

BI-REAGENT DETERMINATION	SENSITIVITY	REFERENCE RANGE	RANGE UP TO
RHEUMA			
ANTISTREPTOLYSIN O	10 IU/mL	250 IU/mL	900 IU/mL
C-REACTIVE PROTEIN	0.4 mg/L	6-8 mg/L	100 mg/L
RHEUMATOID FACTOR	4.2 IU/mL	20 IU/mL	140 IU/mL
INFLAMMATORY STATUS			
ULTRASENSITIVE CRP	0.18 mg/L	6-8 mg/L	50 mg/L
CARDIAC RISK			
CRP CARDIAC MARKER	0.05 mg/L	1-3 mg/L	12.5 mg/L
LIPOPROTEIN (A)	5.0 mg/L	300 mg/L	1200 mg/L
MALNUTRITION STATUS			
RETINOL-BINDING PROTEIN	1.0 mg/L	30-60 mg/L	80 mg/L
RENAL PROFILE			
CYSTATIN C	0.05 mg/L	0.59-1.03 mg/L	10 mg/L
α ₁ -MICROGLOBULIN	1.0 mg/L	10 mg/L	90 mg/L
β ₂ -MICROGLOBULIN	0.2 mg/L	0.8-2.4 mg/L	12 mg/L
IRON METABOLISM			
FERRITIN	5.2 ng/mL	30-300 ng/mL	500 ng/mL
ALLERGY			
IMMUNOGLOBULIN E	10 IU/mL	0-200 IU/mL	1500 IU/mL
DIABETIC PROFILE			
MICROALBUMIN	3.0 mg/L	0-20 µg/min	250 mg/L
MONDREAGENT DETERMINATION	SENSITIVITY	REFERENCE RANGE	LINEAR RANGE
RHEUMA			
ANTISTREPTOLYSIN O	15 IU/mL	250 IU/mL	400 IU/mL
C-REACTIVE PROTEIN	2.0 mg/L	6-8 mg/L	90 mg/L
RHEUMATOID FACTOR	10 IU/mL	20 IU/mL	80 IU/mL
CARDIAC RISK			
CRP CARDIAC MARKER	0.05 mg/L	1-3 mg/L	6.5 mg/L
LIPOPROTEIN (A)	15 mg/L	300 mg/L	400 mg/L
RENAL PROFILE			
β ₂ -MICROGLOBULIN	0.2 mg/L	0.8-2.4 mg/L	15 mg/L
IRON METABOLISM			
FERRITIN	5.0 ng/mL	30-300 ng/mL	300 ng/mL
DIABETIC PROFILE			
MICROALBUMIN	5.0 mg/L	0-20 µg/min	125 mg/L


PROZONE >	ASSAY TIME	CALIBRATION	CONTROLS	R ₁ + R ₂	PAGE
1 500 IU/mL	7 min	Multipoint	2	6L + 1L	6
200 mg/L	5 min	Multipoint	2	5L + 1L	7
1 500 IU/mL	7 min	Multipoint	2	7L + 1L	8
120 mg/L	5 min	Multipoint	2	5L + 1L	9
900 mg/L	5 min	Multipoint	1	5L + 1L	10
2 250 mg/L	5 min	Multipoint	1	6L + 1L	11
160 mg/L	5 min	Multipoint	2	5L + 1L	12
16 mg/L	5 min	Multipoint	2	5L + 1L	13
210 mg/L	7 min	Multipoint	1	5L + 1L	14
100 mg/L	4 min	Multipoint	2	5L + 1L	15
5 500 ng/mL	4 min	Multipoint	2	3L + 1L	16
12 000 IU/mL	4 min	Multipoint	2	3L + 1L	17
500 mg/L	6 min	Multipoint	2	5L + 1L	18
PROZONE >	ASSAY TIME	CALIBRATION	CONTROLS	R ₁ + R ₂	PAGE
1 500 IU/mL	2 min	1 point	1	9L + 1L	19
430 mg/L	2 min	1 point	1	9L + 1L	20
1 500 IU/mL	2.5 min	1 point	1	9L + 1L	21
45 mg/L	5 min	1 point	1	9L + 1L	22
2 200 mg/L	4 min	1 point	1	7L + 1L	23
100 mg/L	5 min	1 point	1	5L + 1L	24
7 000 ng/mL	5 min	1 point	1	2.5L + 1L	25
500 mg/L	4 min	1 point	1	9L + 1L	26

- Range of analysis up to 900 IU/mL
- Excellent precisions (n=10):

	1	2	3
Concentration IU/mL	85	200	435
Intra-assay CV %	3.8	2.8	2.7
Inter-assay CV %	4.6	2.7	2.1

- Sensitivity: least detectable dose (n=20) < 10 IU/mL
- Prozone > 1500 IU/mL

GENERAL APPLICATION



Volume, time and wavelength are recommended and may be adjusted depending on the analyzer features.

ANTISTREPTOLYSIN O

MAIN FEATURES

- Latex enhanced immunoturbidimetric method.
- Easy to program on any type of automated analyzers.
- Correlation coefficient of 0.99 with another commercial procedure.
- Long self-life up to 2 years, liquid reagents ready to use.
- Traceable to the International Standard for Antistreptolysin-O (WHO).

DIAGNOSTIC RELEVANCE

- Antistreptolysin O (ASO) testing is used for the diagnosis of non-suppurative complications of infections caused by streptococcal pathogens: acute rheumatic fever or acute poststreptococcal glomerulonephritis.

MEASURING METHOD

- ASO and latex particles bound streptolysin O reaction.
- ASO values determined turbidimetrically.


BI-REAGENT

- Range of analysis up to 100 mg/L
- Excellent precisions (n=10):

	1	2	3
Concentration mg/L	5.5	35	60
Intra-assay CV %	2.6	1.9	1.7
Inter-assay CV %	5.0	3.2	5.5

- Sensitivity: least detectable dose (n=20) < 0.4 mg/L
- Prozone > 200 mg/L

GENERAL APPLICATION



Volume, time and wavelength are recommended and may be adjusted depending on the analyzer features.

C-REACTIVE PROTEIN

MAIN FEATURES

- Latex enhanced immunoturbidimetric method.
- Easy to program on any type of automated analyzers.
- Correlation coefficient of 0.99 with another commercial procedure.
- Long self-life up to 2 years, liquid reagents ready to use.
- Traceable to CRM 470 International Standard.

DIAGNOSTIC RELEVANCE

- Elevated C-reactive Protein (CRP) has been demonstrated in nearly all bacterial and fungal infections.
- In addition, it has been shown to be increased in other diseases as neoplasia, and rheumatic diseases as well as in major surgery.

MEASURING METHOD

- CRP in the sample reaction with latex-covalently bound antibodies against human CRP.
- CRP values determined turbidimetrically.


biol latex
WWW.BIOLATEX.COM

- Range of analysis up to 140 IU/mL
- Excellent precisions (n=10):

	1	2	3
Concentration IU/mL	40	70	125
Intra-assay CV %	2.3	1.6	2.3
Inter-assay CV %	5.6	2.3	

- Sensitivity: least detectable dose (n=20) < 4.2 IU/mL
- Prozone > 1500 IU/mL

GENERAL APPLICATION



Volume, time and wavelength are recommended and may be adjusted depending on the analyzer features.

RHEUMATOID FACTOR

MAIN FEATURES

- Latex enhanced immunoturbidimetric method.
- Easy to program on any type of automated analyzers.
- Correlation coefficient of 0.97 with another commercial procedure.
- Long self-life up to 2 years, liquid reagents ready to use.
- Standardized against the International Reference Preparation of Rheumatoid Arthritis Serum (WHO).

DIAGNOSTIC RELEVANCE

- Rheumatoid Factors (RF) determinations are clinically important for the diagnosis, prognosis and assessment of therapeutic efficacy of rheumatoid arthritis.

MEASURING METHOD


- IgM-anti-IgG (RF) and latex-covalently bound human IgG reaction.
- RF values determined turbidimetrically.

- Range of analysis up to 50 mg/L
- Excellent precisions (n=10):

	1	2	3
Concentration mg/L	1.75	5.5	40
Intra-assay CV %	5.2	3.2	2.3
Inter-assay CV %	4.2	2.0	5.6

- Sensitivity: least detectable dose (n=20) < 0.18 mg/L
- Prozone > 120 mg/L

GENERAL APPLICATION



Volume, time and wavelength are recommended and may be adjusted depending on the analyzer features.

ULTRASENSITIVE CRP

MAIN FEATURES

- Latex enhanced immunoturbidimetric method.
- Easy to program on any type of automated analyzers.
- Correlation coefficient of 1.00 with another commercial procedure.
- Long self-life up to 2 years, liquid reagents ready to use.
- Traceable to CRM 470 International Standard.

DIAGNOSTIC RELEVANCE

- Elevated C-reactive Protein (CRP) has been demonstrated in nearly all bacterial and fungal infections.
- In addition, it has been shown to be increased in other diseases as neoplasia, and rheumatic diseases as well as in major surgery.

MEASURING METHOD


- CRP in the sample reaction with latex-covalently bound antibodies against human CRP.
- CRP values determined turbidimetrically.

- Range of analysis up to 12.5 mg/L
- Excellent precisions (n=10):

	1	2	3
Concentration mg/L	0.6	1.1	2.1
Intra-assay CV %	9.2	2.6	2.8
Inter-assay CV %	6.8	3.4	3.2

- Sensitivity: least detectable dose (n=20) < 0.05 mg/L
- Prozone > 900 mg/L

GENERAL APPLICATION



Volume, time and wavelength are recommended and may be adjusted depending on the analyzer features.

CRP CARDIAC MARKER

MAIN FEATURES

- Latex enhanced immunoturbidimetric method.
- Easy to program on any type of automated analyzers.
- Correlation coefficient of 0.99 with another commercial procedure.
- Long self-life up to 2 years, liquid reagents ready to use.
- Traceable to CRM 470 International Standard.

DIAGNOSTIC RELEVANCE

- C-reactive Protein (CRP) is a marker of inflammation that seems to be a strong predictor of cardiovascular events and has been used to predict incident myocardial infarction, stroke, peripheral arterial disease and sudden cardiac death.
- The addition of CRP to standard cholesterol evaluation may thus provide a simple and inexpensive method to improve global risk prediction and compliance with preventive approaches.

MEASURING METHOD

- CRP in the sample reaction with latex-covalently bound antibodies against human CRP.
- CRP values determined turbidimetrically.


BI-REAGENT

- Range of analysis up to 1200 mg/L
- Excellent precisions (n=10):

	1	2	3
Concentration mg/L	160	479	715
Intra-assay CV %	4.7	4.5	3.4
Inter-assay CV %	5.1	3.9	3.5

- Sensitivity: least detectable dose (n=20) < 5.0 mg/L
- Prozone > 2250 mg/L

GENERAL APPLICATION



Volume, time and wavelength are recommended and may be adjusted depending on the analyzer features.

LIPOPROTEIN (A)

MAIN FEATURES

- Latex enhanced immunoturbidimetric method.
- Easy to program on any type of automated analyzers.
- Correlation coefficient of 0.96 with another commercial procedure.
- Long self-life up to 2 years, liquid reagents ready to use.
- Standardized against WHO reference material SRM2B.

DIAGNOSTIC RELEVANCE

- High lipoprotein(a) [Lp(a)] concentration is associated with risk for cardiovascular disease and also associated with an increased risk of coronary heart disease.
- Quantitation of Lp(a) is important for identification of individuals at risk for developing atherosclerosis.

MEASURING METHOD

- Lp(a) in the sample and latex-covalently bound antibodies against human Lp(a) reaction.
- Lp(a) values determined turbidimetrically.

biolatrix
WWW.BIOLATEX.COM


BI-REAGENT

- Range of analysis up to 80 mg/L
- Excellent precisions (n=10):

	1	2	3
Concentration mg/L	17	33	60
Intra-assay CV %	2.2	0.8	1.5
Inter-assay CV %	2.7	2.1	3.7

- Sensitivity: least detectable dose (n=20) < 1.0 mg/L
- Prozone > 160 mg/L

GENERAL APPLICATION



Volume, time and wavelength are recommended and may be adjusted depending on the analyzer features.

RETINOL-BINDING PROTEIN

MAIN FEATURES

- Latex enhanced immunoturbidimetric method.
- Easy to program on any type of automated analyzers.
- Correlation coefficient of 0.99 with another commercial procedure.
- Long self-life up to 2 years, liquid reagents ready to use.
- Standardized against a highly purified material.

DIAGNOSTIC RELEVANCE

- Retinol-binding Protein (RBP) concentrations are markedly diminished in malnutrition and other conditions.
- Due to its short half live may be suitable for monitoring the nutritional status and efficacy of parenteral nutrition.

MEASURING METHOD

- RBP and latex-covalently bound antibodies against human RBP reaction.
- RBP values determined turbidimetrically.


bio:latex
WWW.BIOLATEX.COM

- Range of analysis up to 10 mg/L
- Excellent precisions (n=10):

	1	2	3	4
Concentration mg/L	0.85	1.5	3.0	5.0
Intra-assay CV %	1.0	1.2	1.0	0.7
Inter-assay CV %	3.4	2.2	1.5	2.9

- Sensitivity: least detectable dose (n=20) < 0.05 mg/L
- Prozone > 16 mg/L

GENERAL APPLICATION



Volume, time and wavelength are recommended and may be adjusted depending on the analyzer features.

CYSTATIN C

MAIN FEATURES

- Latex enhanced immunoturbidimetric method.
- Easy to program on any type of automated analyzers.
- Correlation coefficient of 0.99 with another commercial procedure.
- Long self-life up to 2 years, liquid reagents ready to use.
- Standardized against IFCC human serum ERM-DA471.

DIAGNOSTIC RELEVANCE

- Cystatin C is an excellent indicator of GFR (glomerular filtration rate) as its concentration is almost exclusively determined by the GFR.
- Cystatin C is also acknowledged as a marker of elevated risk of death from cardiovascular complications – myocardial infarction and stroke.

MEASURING METHOD

- Cystatin C and latex-covalently bound antibodies against human Cystatin C reaction.
- Cystatin C values determined turbidimetrically.


BI-REAGENT

- Range of analysis up to 90 mg/L
- Excellent precisions (n=10):

	1	2	3
Concentration mg/L	16	28	40
Intra-assay CV %	2.4	3.9	3.1
Inter-assay CV %	2.8	4.3	5.2

- Sensitivity: least detectable dose (n=20) < 1.0 mg/L
- Prozone > 210 mg/L

GENERAL APPLICATION



Volume, time and wavelength are recommended and may be adjusted depending on the analyzer features.

α_1 -MICROGLOBULIN

MAIN FEATURES

- Latex enhanced immunoturbidimetric method.
- Easy to program on any type of automated analyzers.
- Correlation coefficient of 0.98 with another commercial procedure.
- Long self-life up to 2 years, liquid reagents ready to use.
- Standardized with reference to highly purified proteins preparation.

DIAGNOSTIC RELEVANCE

- Elevated concentrations of low molecular weight proteins in urine such as α_1 -microglobulin indicates tubular damage, which can occur after heavy metal exposure or in the course of advanced diabetic nephropathy, nephritis or others pathologies.

MEASURING METHOD

- Antibody to human α_1 -microglobulin bound to latex particles reaction with α_1 -microglobulin in a sample.
- α_1 -microglobulin values determined turbidimetrically.


● ● ● BI-REAGENT

- Range of analysis up to 12 mg/L
- Excellent precisions (n=10):

	1	2	3
Concentration mg/L	1.8	5.8	11.5
Intra-assay CV %	2.3	2.2	2.1
Inter-assay CV %	4.1	3.9	7.5

- Sensitivity: least detectable dose (n=20) < 0.2 mg/L
- Prozone > 100 mg/L

GENERAL APPLICATION



Volume, time and wavelength are recommended and may be adjusted depending on the analyzer features.

β₂-MICROGLOBULIN

MAIN FEATURES

- Latex enhanced immunoturbidimetric method.
- Easy to program on any type of automated analyzers.
- Correlation coefficient of 0.99 with another commercial procedure.
- Long self-life up to 2 years, liquid reagents ready to use.
- Standardized with reference to highly purified proteins preparation.

DIAGNOSTIC RELEVANCE

- Increased urinary excretion of β₂-microglobulin is an indicator of renal tubular disorders used to detect early nephrotoxicity in patients treated with gentamicin and other nephrotoxic drugs.
- Elevated serum levels has been reported as a useful marker of Acquired Immune Deficiency Syndrome, in myeloma patients and in a variety of diseases including carcinomas and lymphoid tumours and inflammatory and autoimmune diseases.

MEASURING METHOD

- β₂-microglobulin in the sample reaction with latex-covalently bound antihuman β₂-microglobulin antibodies.
- β₂-microglobulin values determined turbidimetrically.


● ● ● *biol latex*
WWW.BIOLATEX.COM

- Range of analysis up to 500 ng/mL
- Excellent precisions (n=10):

	1	2
Concentration ng/mL	180	285
Intra-assay CV %	1.6	1.6
Inter-assay CV %	3.3	3.2

- Sensitivity: least detectable dose (n=20) < 5.2 ng/mL
- Prozone > 5500 ng/mL

GENERAL APPLICATION



Volume, time and wavelength are recommended and may be adjusted depending on the analyzer features.

FERRITIN

MAIN FEATURES

- Latex enhanced immunoturbidimetric method.
- Easy to program on any type of automated analyzers.
- Correlation coefficient of 1.00 with another commercial procedure.
- Long self-life up to 18 months, liquid reagents ready to use.
- Traceable to WHO 80/578 International Standard.

DIAGNOSTIC RELEVANCE

- Ferritin concentration is correlated with the quantity of available iron stored in the body so its determination is used for diagnosis and monitoring of iron deficiency and iron overload.
- Concentrations of Ferritin are found to be elevated in patients with infections, inflammation or in hepatic or chronic renal diseases.

MEASURING METHOD


- Ferritin in the sample and latex-covalently bound antibodies against human Ferritin reaction.
- Ferritin values determined turbidimetrically.

- Range of analysis up to 1500 IU/mL
- Excellent precisions (n=10):

	1	2	3
Concentration IU/mL	40	80	140
Intra-assay CV %	3.3	1.6	1.3
Inter-assay CV %	4.7	3.4	1.5

- Sensitivity: least detectable dose (n=20) < 10 IU/mL
- Prozone > 12000 IU/mL

GENERAL APPLICATION



Volume, time and wavelength are recommended and may be adjusted depending on the analyzer features.

IMMUNOGLOBULIN E

MAIN FEATURES

- Latex enhanced immunoturbidimetric method.
- Easy to program on any type of automated analyzers.
- Correlation coefficient of 1.00 with another commercial procedure.
- Long self-life up to 1 year, liquid reagents ready to use.
- Standardized against IRP 75/502.

DIAGNOSTIC RELEVANCE

- Immunoglobulin (IgE) determinations are indicated in the diagnosis and monitoring of allergic diseases.
- Elevated IgE levels also occur in parasitosis and immunodeficiency syndromes.
- In infants and small children with recurrent respiratory tract diseases the determination of IgE is of prognostic relevance, and also in some myelomas of IgE type.

MEASURING METHOD

- Anti-IgE antibodies covalently bound to latex particles reaction with the IgE in the sample.
- IgE values determined turbidimetrically.


BI-REAGENT

- Range of analysis up to 250 mg/L
- Excellent precisions (n=10):

	1	2	3
Concentration mg/L	20	55	115
Intra-assay CV %	0.8	0.6	1.1
Inter-assay CV %	6.3	1.9	4.1

- Sensitivity: least detectable dose (n=20) < 3.0 mg/L
- Prozone > 500 mg/L

GENERAL APPLICATION



Volume, time and wavelength are recommended and may be adjusted depending on the analyzer features.

MICROALBUMIN

MAIN FEATURES

- Latex enhanced immunoturbidimetric method.
- Easy to program on any type of automated analyzers.
- Correlation coefficient of 0.98 with another commercial procedure.
- Long self-life up to 2 years, liquid reagents ready to use.
- Traceable to CRM 470 International Standard.

DIAGNOSTIC RELEVANCE

- Elevations in the albumin excretion has been used as a predictor of nephropathy and cardiovascular disease in diabetic patients.
- It has also been associated with hypertension and increased risk of cardiovascular disease in non-diabetic patients.

MEASURING METHOD

- Sample albumin reaction with latex enhanced albumin antibody.
- Microalbumin concentration determined turbidimetrically.


biolatrix
WWW.BIOLATEX.COM

- Range of analysis up to 900 IU/mL (linear up to 400 IU/mL)
- Excellent precisions (n=10):

	1	2	3
Concentration IU/mL	100	200	430
Intra-assay CV %	2.9	3.1	2.9
Inter-assay CV %	2.9	3.6	4,4

- Sensitivity: least detectable dose (n=20) < 15 IU/mL
- Prozone > 1500 IU/mL

GENERAL APPLICATION



(*) Working Reagent (WR)
1 part of reagent latex +
9 parts of reagent buffer

ANTISTREPTOLYSIN O

MAIN FEATURES

- Latex enhanced immunoturbidimetric monoreagent method.
- Easy to program on any type of automated analyzers.
- Correlation coefficient of 0.99 with another commercial procedure.
- Long self-life up to 2 years, liquid reagents ready to use.
- Traceable to the International Standard for Antistreptolysin-O (WHO).

DIAGNOSTIC RELEVANCE

- Antistreptolysin O (ASO) testing is used for the diagnosis of non-suppurative complications of infections caused by streptococcal pathogens: acute rheumatic fever or acute poststreptococcal glomerulonephritis.

MEASURING METHOD


- ASO and latex particles bound streptolysin O reaction.
- ASO values determined photometrically.

- Range of analysis up to 90 mg/L (linear up to 90 mg/L)
- Excellent precisions (n=10):

	1	2	3
Concentration mg/L	8	18	33
Intra-assay CV %	2.8	3.2	4.2
Inter-assay CV %	3.9	4.1	4.6

- Sensitivity: least detectable dose (n=20) < 2.0 mg/L
- Prozone > 430 mg/L

GENERAL APPLICATION



(*) Working Reagent (WR)
1 part of reagent latex +
9 parts of reagent buffer

C-REACTIVE PROTEIN

MAIN FEATURES

- Latex enhanced immunoturbidimetric monoreagent method.
- Easy to program on any type of automated analyzers.
- Correlation coefficient of 0.99 with another commercial procedure.
- Long self-life up to 2 years, liquid reagents ready to use.
- Traceable to CRM 470 International Standard.

DIAGNOSTIC RELEVANCE

- Elevated C-reactive Protein (CRP) has been demonstrated in nearly all bacterial and fungal infections.
- In addition, it has been shown to be increased in other diseases as neoplasia, and rheumatic diseases as well as in major surgery.

MEASURING METHOD


- CRP in the sample reaction with latex-covalently bound antibodies against human CRP.
- CRP values determined photometrically.

- Range of analysis up to 200 IU/mL (linear up to 80 IU/mL)
- Excellent precisions (n=10):

	1	2	3
Concentration IU/mL	17	70	120
Intra-assay CV %	4.7	2.6	2.1
Inter-assay CV %	2.5	3.8	2.3

- Sensitivity: least detectable dose (n=20) < 10 IU/mL
- Prozone > 1500 IU/mL

GENERAL APPLICATION



(*) Working Reagent (WR)
1 part of reagent latex +
9 parts of reagent buffer

RHEUMATOID FACTOR

MAIN FEATURES

- Latex enhanced immunoturbidimetric monoreagent method.
- Easy to program on any type of automated analyzers.
- Long self-life up to 2 years, liquid reagents ready to use.
- Standardized against the International Reference Preparation of Rheumatoid Arthritis Serum (WHO).

DIAGNOSTIC RELEVANCE

- Rheumatoid Factors (RF) determinations are clinically important for the diagnosis, prognosis and assessment of therapeutic efficacy of rheumatoid arthritis.

MEASURING METHOD


- IgM-anti-IgG (RF) and latex-covalently bound human IgG reaction.
- RF values determined photometrically.

- Range of analysis up to 12.5 mg/L (linear up to 6.5 mg/L)
- Excellent precisions (n=10):

	1	2	3
Concentration mg/L	0.75	1.5	3.0
Intra-assay CV %	2.3	2.4	1.6
Inter-assay CV %	2.5	2.0	1.6

- Sensitivity: least detectable dose (n=20) < 0.05 mg/L
- Prozone > 45 mg/L

GENERAL APPLICATION



(*) Working Reagent (WR)
1 part of reagent latex +
9 parts of reagent buffer

CRP CARDIAC MARKER

MAIN FEATURES

- Latex enhanced immunoturbidimetric monoreagent method.
- Easy to program on any type of automated analyzers.
- Correlation coefficient of 0.98 with another commercial procedure.
- Long self-life up to 2 years, liquid reagents ready to use.
- Traceable to CRM 470 International Standard.

DIAGNOSTIC RELEVANCE

- C-reactive Protein (CRP) is a marker of inflammation that seems to be a strong predictor of cardiovascular events and has been used to predict incident myocardial infarction, stroke, peripheral arterial disease and sudden cardiac death.
- The addition of CRP to standard cholesterol evaluation may thus provide a simple and inexpensive method to improve global risk prediction and compliance with preventive approaches.

MEASURING METHOD


- CRP in the sample reaction with latex-covalently bound antibodies against human CRP.
- CRP values determined photometrically.

- Range of analysis up to 900 mg/L (linear up to 400 mg/L)
- Excellent precisions (n=10):

	1	2	3
Concentration mg/L	90	170	325
Intra-assay CV %	1.7	2.0	4.4
Inter-assay CV %	3.3	1.7	4.8

- Sensitivity: least detectable dose (n=20) < 15 mg/L
- Prozone > 2200 mg/L

GENERAL APPLICATION



(*) Working Reagent (WR)
1 part of reagent latex +
7 parts of reagent buffer

LIPOPROTEIN (A)

MAIN FEATURES

- Latex enhanced immunoturbidimetric monoreagent method.
- Easy to program on any type of automated analyzers.
- Correlation coefficient of 0.99 with another commercial procedure.
- Long self-life up to 2 years, liquid reagents ready to use.
- Standardized against WHO reference material SRM2B.

DIAGNOSTIC RELEVANCE

- High lipoprotein(a) [Lp(a)] concentration is associated with risk for cardiovascular disease and also associated with an increased risk of coronary heart disease.
- Quantitation of Lp(a) is important for identification of individuals at risk for developing atherosclerosis.

MEASURING METHOD


- Lp(a) in the sample and latex-covalently bound antibodies against human Lp(a) reaction.
- Lp(a) values determined photometrically.

- Range of analysis up to 20 mg/L (linear up to 15 mg/L)
- Excellent precisions (n=10):

	1	2	3
Concentration mg/L	1.75	8	15
Intra-assay CV %	3.2	1.3	1.9
Inter-assay CV %	3.0	1.8	4.2

- Sensitivity: least detectable dose (n=20) < 0.2 mg/L
- Prozone > 100 mg/L

GENERAL APPLICATION



(*) Working Reagent (WR)
1 part of reagent latex +
5 parts of reagent buffer

β₂-MICROGLOBULIN

MAIN FEATURES

- Latex enhanced immunoturbidimetric monoreagent method.
- Easy to program on any type of automated analyzers.
- Correlation coefficient of 1.00 with another commercial procedure.
- Long self-life up to 2 years, liquid reagents ready to use.
- Standardized with reference to highly purified proteins preparation.

DIAGNOSTIC RELEVANCE

- Increased urinary excretion of β₂-microglobulin is an indicator of renal tubular disorders used to detect early nephrotoxicity in patients treated with gentamicin and other nephrotoxic drugs.
- Elevated serum levels has been reported as a useful marker of Acquired Immune Deficiency Syndrome, in myeloma patients and in a variety of diseases including carcinomas and lymphoid tumours and inflammatory and autoimmune diseases.

MEASURING METHOD


- β₂-microglobulin in the sample reaction with latex-covalently bound antihuman β₂-microglobulin antibodies.
- β₂-microglobulin values determined photometrically.

- Range of analysis up to 500 ng/mL (linear up to 300 ng/mL)
- Excellent precisions (n=10):

	1	2	3
Concentration ng/mL	37	170	270
Intra-assay CV %	3.4	1.4	1.0
Inter-assay CV %	3.9	1.7	1.3

- Sensitivity: least detectable dose (n=20) < 5.0 ng/mL
- Prozone > 7000 ng/mL

GENERAL APPLICATION



(*) Working Reagent (WR)
1 part of reagent latex +
2.5 parts of reagent buffer

FERRITIN

MAIN FEATURES

- Latex enhanced immunoturbidimetric monoreagent method.
- Easy to program on any type of automated analyzers.
- Correlation coefficient of 0.95 with another commercial procedure.
- Long self-life up to 18 months, liquid reagents ready to use.
- Traceable to WHO 80/578 International Standard.

DIAGNOSTIC RELEVANCE

- Ferritin concentration is correlated with the quantity of available iron stored in the body so its determination is used for diagnosis and monitoring of iron deficiency and iron overload.
- Concentrations of Ferritin are found to be elevated in patients with infections, inflammation or in hepatic or chronic renal diseases.

MEASURING METHOD


- Ferritin in the sample and latex-covalently bound antibodies against human Ferritin reaction.
- Ferritin values determined photometrically.

- Range of analysis up to 250 mg/L (linear up to 125 mg/L)
- Excellent precisions (n=10):

	1	2	3
Concentration mg/L	30	60	160
Intra-assay CV %	3.4	3.7	1.9
Inter-assay CV %	4.5	2.4	4.5

- Sensitivity: least detectable dose (n=20) < 5.0 mg/L
- Prozone > 500 mg/L

GENERAL APPLICATION



(*) Working Reagent (WR)
1 part of reagent latex +
9 parts of reagent buffer

MICROALBUMIN

MAIN FEATURES

- Latex enhanced immunoturbidimetric monoreagent method.
- Easy to program on any type of automated analyzers.
- Correlation coefficient of 0.98 with another commercial procedure.
- Long self-life up to 2 years, liquid reagents ready to use.
- Traceable to CRM 470 International Standard.

DIAGNOSTIC RELEVANCE

- Elevations in the albumin excretion has been used as a predictor of nephropathy and cardiovascular disease in diabetic patients.
- It has also been associated with hypertension and increased risk of cardiovascular disease in non-diabetic patients.


MEASURING METHOD

- Sample albumin reaction with latex enhanced albumin antibody.
- Microalbumin concentration determined photometrically.



⋮ *biolatrix*

WWW.BIOLATEX.COM



BIOLATEX, S.L.

POL. INDUSTRIAL LA PORTALADA
CALAHORRA, 4 · PAB. 6
26006 LOGROÑO (SPAIN)
TEL.: (34) 941 270 047
FAX: (34) 941 271 077
INFO@BIOLATEX.COM

